




# ***C-VIEW ENDOSCOPE***


InterScience, Inc. is applying its C-View optical technology to medical endoscopes for minimally invasive surgery applications. The C-View optical technology is capable of providing an ultrawide field of view of approximately  $\pm 120^\circ$  about the optical axis. The technology under development can be broadly applied to rigid, semi-rigid, flexible, distal sensor and capsule endoscopy applications.

The ultrawide field of view would have significant impact in medical fields such as:

- Otolaryngology
- Gastroenterology
- Urology
- Gynecology
- Minimally Invasive Surgery (Laparoscopy)



**C-View FOV**



**Standard FOV**